





COL1A1 Monoclonal Antibody

Product Code	CSB-MA955908
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02452
Immunogen	Synthetic Peptide of Collagen I
Raised In	Mouse
Species Reactivity	Human, Rat
Tested Applications	ELISA, IHC; Recommended dilution: IHC: 1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	COL1A1


Image



Immunohistochemical analysis of paraffin-embedded Human skin. COL1A1 Antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Rat-kidney tissue. COL1A1 Antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human-lung-cancer tissue. COL1A1 Antibody was diluted at 1:200.