





COL1A1 Monoclonal Antibody

Product Code	CSB-MA955908
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02452
Immunogen	Synthetic Peptide of Collagen I
Raised In	Mouse
Species Reactivity	Human, Rat
Tested Applications	ELISA, IHC; Recommended dilution: IHC: 1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COL1A1
Clone No.	6B7


Image



Immunohistochemical analysis of paraffin-embedded Human skin. COL1A1 Antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Rat-kidney tissue. COL1A1 Antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human-lung-cancer tissue. COL1A1 Antibody was diluted at 1:200.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.