



Glu-Glu Tag Monoclonal Antibody

Product Code	CSB-MA418664
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Research Use	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	N/A
Immunogen	Synthetic Peptide of Glu-Glu Tag
Raised In	Mouse
Species Reactivity	All
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Product Type	Monoclonal Antibody
Clone No.	24G11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.