





TTR Monoclonal Antibody

Product Code	CSB-MA111203
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Research Use	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02766
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB: 1:500-1:5000, IHC: 1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Alias	Transthyretin (ATTR) (Prealbumin) (TBPA)
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	27H2


Image



Immunohistochemical analysis of paraffin-embedded Human Liver Carcinoma Tissue using TTR Mouse mAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Thyroid Tissue using TTR Mouse mAb diluted at 1:200



Western blot analysis of Human Serum using TTR Mouse mAb diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.