





# TTR Monoclonal Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-MA362195   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Research Use</b>        | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.                  |
| <b>Uniprot No.</b>         | P02766   |
| <b>Immunogen</b>           | Recombinant Protein  |
| <b>Raised In</b>           | Mouse  |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA, WB, IHC; Recommended dilution: WB: 1:500-1:5000, IHC: 1:50-1:500                                    |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.                               |
| <b>Purification Method</b> | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| <b>Alias</b>               | Transthyretin (ATTR) (Prealbumin) (TBPA)   |
| <b>Product Type</b>        | Monoclonal Antibody  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |

## Image



Immunohistochemical analysis of paraffin-embedded Human Thyroid Tissue using TTR Mouse mAb diluted at 1:200



Western blot analysis of Human Serum using TTR Mouse mAb diluted at 1:2000