





COL4A1 Monoclonal Antibody

| | |
|----------------------------|--|
| Product Code | CSB-MA457742 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Research Use | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | p02462 |
| Immunogen | Synthetic Peptide |
| Raised In | Mouse |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,IHC;Recommended dilution:IHC:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| Alias | Collagen alpha-1(IV) chain [Cleaved into: Arresten] |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Clone No. | 7D5 |

Image



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using Collagen IV Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Skin Tissue using Collagen IV Mouse mAb diluted at 1:200.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.