



# Goat Anti-Mouse IgG(H+L) Antibody;HRP conjugated

<b>Product Code</b>	CSB-PA644737
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Goat
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA,WB,IHC,IF;Recommended dilution: WB:1:5000-20000,IHC:1:5000-10000,IF:1:5000-10000,ELISA :1:10000-20000
<b>Form</b>	Liquid
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
<b>Product Type</b>	Secondary Antibody
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.