



Goat Anti-Rabbit IgG(H+L) Antibody;HRP conjugated

Product Code	CSB-PA489724
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Rabbit
Tested Applications	ELISA,WB,IHC,IF;Recommended dilution: WB:1:5000-20000,IHC:1:5000-10000,IF:1:5000-10000,ELISA :1:10000-20000
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
Product Type	Secondary Antibody
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.