



# Goat Anti-Rabbit IgG(H+L) Antibody;HRP conjugated

<b>Product Code</b>	CSB-PA489724
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Goat
<b>Species Reactivity</b>	Rabbit
<b>Tested Applications</b>	ELISA, WB, IHC, IF; Recommended dilution: WB: 1:5000-20000, IHC: 1:5000-10000, IF: 1:5000-10000, ELISA : 1:10000-20000
<b>Form</b>	Liquid
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
<b>Product Type</b>	Secondary Antibody
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)