



Goat Anti-Mouse IgG(H+L) Antibody;FITC conjugated

Product Code	CSB-PA980279
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Mouse
Tested Applications	IHC,IF,FC;Recommended dilution:IHC:1:5000-10000,IF:1:5000-1:10000,FC:50-200
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
Product Type	Secondary Antibody
Immunogen Species	Mus musculus (Mouse)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.