



Goat Anti-Rabbit IgG(H+L) Antibody;FITC conjugated

Product Code	CSB-PA340892
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Rabbit
Tested Applications	IHC,IF,FC; IHC:1/5000-1/10000, IF:1/5000-1/10000, FC:50-200
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
Product Type	Secondary Antibody
Immunogen Species	Oryctolagus cuniculus (Rabbit)