



Goat Anti-Rabbit IgG(H+L) Antibody;FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA340892 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Raised In | Goat |
| Species Reactivity | Rabbit |
| Tested Applications | IHC,IF,FC; IHC:1/5000-1/10000, IF:1/5000-1/10000, FC:50-200 |
| Form | Liquid |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen. |
| Product Type | Secondary Antibody |
| Immunogen Species | Oryctolagus cuniculus (Rabbit) |