



Goat Anti-Rabbit IgG(H+L) Antibody;Biotin conjugated

Product Code	CSB-PA153887
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Rabbit
Tested Applications	ELISA, WB, IHC, IF, IP, FC; Recommended dilution: WB: 1:5000-20000, IHC: 1:200-2000, IF: 1:200-2000, IP: 1:200-1:2000, FCM: 1:200-2000, ELISA: 1:5000-20000
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
Product Type	Secondary Antibody
Immunogen Species	Oryctolagus cuniculus (Rabbit)