



Goat Anti-Rabbit IgG(H+L) Antibody;Biotin conjugated

Product Code	CSB-PA153887
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Rabbit
Tested Applications	ELISA,WB,IHC,IF,IP,FC;Recommended dilution:WB:1:5000-20000,IHC:1:200-2000,IF:1:200-2000,IP:1:200-1:2000,FCM:1:200-2000,ELISA:1:5000-20000
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
Product Type	Secondary Antibody
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.