



Goat Anti-Mouse IgG(H+L) Antibody

Product Code	CSB-PA477109
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB: 1:5000-20000, IHC: 1:5000-10000, IF: 1:5000-10000, ELISA: 1:10000-20000
Form	Liquid
Storage Buffer	Liquid in PBS
Purification Method	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
Product Type	Secondary Antibody
Immunogen Species	Mus musculus (Mouse)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.