



ABCC4 Antibody, HRP conjugated

Product Code	CSB-PA001063LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15439
Immunogen	Recombinant Human Multidrug resistance-associated protein 4 protein (562-697AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multidrug resistance-associated protein 4 (ATP-binding cassette sub-family C member 4) (MRP/cMOAT-related ABC transporter) (Multi-specific organic anion transporter B) (MOAT-B), ABCC4, MRP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ABCC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.