





ABCC9 Antibody

Product Code	CSB-PA001067LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60706
Immunogen	Recombinant Human ATP-binding cassette sub-family C member 9 protein (560-669AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family C member 9 (Sulfonylurea receptor 2), ABCC9, SUR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ABCC9

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA001067LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA001067LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)