



AHRR Antibody, FITC conjugated

Product Code	CSB-PA001482LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A9YTQ3
Immunogen	Recombinant Human Aryl hydrocarbon receptor repressor protein (255-506AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aryl hydrocarbon receptor repressor (AhR repressor) (AhRR) (Class E basic helix-loop-helix protein 77) (bHLHe77), AHRR, BHLHE77 KIAA1234
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	AHRR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.