



ATG4B Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA002290LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y4P1 |
| Immunogen | Recombinant Human Cysteine protease ATG4B protein (350-400AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Cysteine protease ATG4B (EC 3.4.22.-) (AUT-like 1 cysteine endopeptidase) (Autophagin-1) (Autophagy-related cysteine endopeptidase 1) (Autophagy-related protein 4 homolog B) (hAPG4B), ATG4B, APG4B AUTL1 KIAA0943 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cancer |
| Target Names | ATG4B |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |