



ATG4B Antibody, Biotin conjugated

Product Code	CSB-PA002290LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4P1
Immunogen	Recombinant Human Cysteine protease ATG4B protein (350-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cysteine protease ATG4B (EC 3.4.22.-) (AUT-like 1 cysteine endopeptidase) (Autophagin-1) (Autophagy-related cysteine endopeptidase 1) (Autophagy-related protein 4 homolog B) (hAPG4B), ATG4B, APG4B AUTL1 KIAA0943
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	ATG4B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.