



BACE1 Antibody, Biotin conjugated

Product Code	CSB-PA002524LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56817
Immunogen	Recombinant Human Beta-secretase 1 protein (46-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-secretase 1 (EC 3.4.23.46) (Aspartyl protease 2) (ASP2) (Asp 2) (Beta-site amyloid precursor protein cleaving enzyme 1) (Beta-site APP cleaving enzyme 1) (Memapsin-2) (Membrane-associated aspartic protease 2), BACE1, BACE KIAA1149
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	BACE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.