



CALCA Antibody, HRP conjugated

Product Code	CSB-PA004434HB01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41547
Immunogen	Recombinant Dog Calcitonin protein (85-116AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcitonin, CALCA, CALC CCALCI
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Research Area	Others
Target Names	CALCA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.