



CBL Antibody, FITC conjugated

Product Code	CSB-PA004578LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22681
Immunogen	Recombinant Human E3 ubiquitin-protein ligase CBL protein (548-847AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase CBL (EC 2.3.2.27) (Casitas B-lineage lymphoma proto-oncogene) (Proto-oncogene c-Cbl) (RING finger protein 55) (RING-type E3 ubiquitin transferase CBL) (Signal transduction protein CBL), CBL, CBL2 RNF55
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CBL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.