





CDKN2C Antibody

Product Code	CSB-PA005094LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42773
Immunogen	Recombinant Human Cyclin-dependent kinase 4 inhibitor C protein (1-68AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase 4 inhibitor C (Cyclin-dependent kinase 6 inhibitor) (p18-INK4c) (p18-INK6), CDKN2C, CDKN6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDKN2C

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA005094LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA005094LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)