




Cpz Antibody

Product Code	CSB-PA005931LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O54858
Immunogen	Recombinant Rat Carboxypeptidase Z protein (21-652AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carboxypeptidase Z (CPZ) (EC 3.4.17.-), Cpz
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Cpz

Image



Western Blot

Positive WB detected in: HEK293 whole cell lysate, A549 whole cell lysate, U87 whole cell lysate

All lanes: CPZ antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 74 kDa

Observed band size: 74, 59 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.