



Cpz Antibody, Biotin conjugated

Product Code	CSB-PA005931LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O54858
Immunogen	Recombinant Rat Carboxypeptidase Z protein (21-652AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carboxypeptidase Z (CPZ) (EC 3.4.17.-), Cpz
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Cpz