



CTH Antibody, FITC conjugated

Product Code	CSB-PA006160LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32929
Immunogen	Recombinant Human Cystathionine gamma-lyase protein (1-70AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cystathionine gamma-lyase (EC 4.4.1.1) (Cysteine-protein sulfhydrase) (Gamma-cystathionase), CTH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CTH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.