




DGAT1 Antibody

Product Code	CSB-PA006817LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75907
Immunogen	Recombinant Human Diacylglycerol O-acyltransferase 1 protein (1-62AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diacylglycerol O-acyltransferase 1 (EC 2.3.1.20) (ACAT-related gene product 1) (Acyl-CoA retinol O-fatty-acyltransferase) (ARAT) (Retinol O-fatty-acyltransferase) (EC 2.3.1.76) (Diglyceride acyltransferase), DGAT1, AGRP1 DGAT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DGAT1

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA006817LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.