


DKC1 Antibody

Product Code	CSB-PA006919LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60832
Immunogen	Recombinant Human H/ACA ribonucleoprotein complex subunit 4 protein (94-170AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	H/ACA ribonucleoprotein complex subunit 4 (EC 5.4.99.-) (CBF5 homolog) (Dyskerin) (Nopp140-associated protein of 57 kDa) (Nucleolar protein NAP57) (Nucleolar protein family A member 4) (snoRNP protein DKC1), DKC1, NOLA4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DKC1

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate, Rat brain tissue, Rat heart tissue, Mouse brain tissue, Mouse heart tissue, Mouse kidney tissue


All lanes: DKC1 antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 58, 48 kDa

Observed band size: 48 kDa



Immunofluorescent analysis of HeLa cells using CSB-PA006919LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.