



DOK3 Antibody, FITC conjugated

Product Code	CSB-PA007019LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L591
Immunogen	Recombinant Human Docking protein 3 protein (1-146AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Docking protein 3 (Downstream of tyrosine kinase 3), DOK3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DOK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.