



DNASE2 Antibody, Biotin conjugated

Product Code	CSB-PA007053LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00115
Immunogen	Recombinant Human Deoxyribonuclease-2-alpha protein (19-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyribonuclease-2-alpha (EC 3.1.22.1) (Acid DNase) (Deoxyribonuclease II alpha) (DNase II alpha) (Lysosomal DNase II) (R31240_2), DNASE2, DNASE2A DNL2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DNASE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.