





DXO Antibody

Product Code	CSB-PA007116LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O77932
Immunogen	Recombinant Human Decapping and exoribonuclease protein (69-266AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Decapping and exoribonuclease protein (DXO) (EC 3.1.13.-) (EC 3.6.1.-) (Dom-3 homolog Z), DXO, DOM3L DOM3Z NG6
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DXO

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA007116LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA007116LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.