





DPYSL4 Antibody

Product Code	CSB-PA007172LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14531
Immunogen	Recombinant Human Dihydropyrimidinase-related protein 4 protein (280-560AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydropyrimidinase-related protein 4 (DRP-4) (Collapsin response mediator protein 3) (CRMP-3) (UNC33-like phosphoprotein 4) (ULIP-4), DPYSL4, CRMP3 ULIP4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DPYSL4


Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA007172LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA007172LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA007172LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.