



DPYSL4 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA007172LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O14531 |
| Immunogen | Recombinant Human Dihydropyrimidinase-related protein 4 protein (280-560AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Dihydropyrimidinase-related protein 4 (DRP-4) (Collapsin response mediator protein 3) (CRMP-3) (UNC33-like phosphoprotein 4) (ULIP-4), DPYSL4, CRMP3 ULIP4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | DPYSL4 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |