



DPYSL4 Antibody, Biotin conjugated

Product Code	CSB-PA007172LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14531
Immunogen	Recombinant Human Dihydropyrimidinase-related protein 4 protein (280-560AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydropyrimidinase-related protein 4 (DRP-4) (Collapsin response mediator protein 3) (CRMP-3) (UNC33-like phosphoprotein 4) (ULIP-4), DPYSL4, CRMP3 ULIP4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DPYSL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.