



EDNRB Antibody, HRP conjugated

Product Code	CSB-PA007404LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24530
Immunogen	Recombinant Human Endothelin receptor type B protein (27-101AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endothelin receptor type B (ET-B) (ET-BR) (Endothelin receptor non-selective type), EDNRB, ETRB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EDNRB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.