



EIF4B Antibody, Biotin conjugated

Product Code	CSB-PA007554LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23588
Immunogen	Recombinant Human Eukaryotic translation initiation factor 4B protein (362-590AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eukaryotic translation initiation factor 4B (eIF-4B), EIF4B
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EIF4B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.