



EN1 Antibody, Biotin conjugated

Product Code	CSB-PA007659LD01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05916
Immunogen	Recombinant Chicken Homeobox protein engrailed-1 protein (1-333AA)
Raised In	Rabbit
Species Reactivity	Chicken
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein engrailed-1 (Gg-En-1) (Homeobox protein en-1), EN1, EN-1
Immunogen Species	Gallus gallus (Chicken)
Research Area	Neuroscience
Target Names	EN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.