



# EPHA5 Antibody, Biotin conjugated

Product Code	CSB-PA007725LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54756
Immunogen	Recombinant Human Ephrin type-A receptor 5 protein (246-514AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin type-A receptor 5 (EC 2.7.10.1) (Brain-specific kinase) (EPH homology kinase 1) (EHK-1) (EPH-like kinase 7) (EK7) (hEK7), EPHA5, BSK EHK1 HEK7 TYRO4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EPHA5