



ERN1 Antibody, Biotin conjugated

Product Code	CSB-PA007795LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75460
Immunogen	Recombinant Human Serine/threonine-protein kinase/endoribonuclease IRE1 protein (493-561AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase/endoribonuclease IRE1 (Endoplasmic reticulum-to-nucleus signaling 1) (Inositol-requiring protein 1) (hIRE1p) (Ire1-alpha) (IRE1a) [Includes: Serine/threonine-protein kinase (EC 2.7.11.1); Endoribonuclease (EC 3.1.26.-)], ERN1, IRE1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ERN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.