



FAAH Antibody, Biotin conjugated

Product Code	CSB-PA007938LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00519
Immunogen	Recombinant Human Fatty-acid amide hydrolase 1 protein (205-342AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fatty-acid amide hydrolase 1 (EC 3.5.1.99) (Anandamide amidohydrolase 1) (Oleamide hydrolase 1), FAAH, FAAH1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FAAH