



FADS1 Antibody, Biotin conjugated

Product Code	CSB-PA007959LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60427
Immunogen	Recombinant Human Fatty acid desaturase 1 protein (1-121AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fatty acid desaturase 1 (EC 1.14.19.-) (Delta(5) fatty acid desaturase) (D5D) (Delta(5) desaturase) (Delta-5 desaturase), FADS1, FADSD5
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	FADS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.