



# FOXE1 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA008811LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O00358
<b>Immunogen</b>	Recombinant Human Forkhead box protein E1 protein (184-365AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Forkhead box protein E1 (Forkhead box protein E2) (Forkhead-related protein FKHL15) (HFKH4) (HNF-3/fork head-like protein 5) (HFKL5) (Thyroid transcription factor 2) (TTF-2), FOXE1, FKHL15 FOXE2 TITF2 TTF2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	FOXE1