




FOXN1 Antibody

Product Code	CSB-PA008829LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15353
Immunogen	Recombinant Human Forkhead box protein N1 protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein N1 (Winged-helix transcription factor nude), FOXN1, RONU WHN
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FOXN1

Image



Western Blot

Positive WB detected in: Raji whole cell lysate, Mouse lung tissue, Mouse stomach tissue


All lanes: FOXN1 antibody at 2.7µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 69 kDa

Observed band size: 69 kDa



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA008829LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA008829LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.