



FOXN1 Antibody, Biotin conjugated

Product Code	CSB-PA008829LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15353
Immunogen	Recombinant Human Forkhead box protein N1 protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein N1 (Winged-helix transcription factor nude), FOXN1, RONU WHN
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FOXN1