



FURIN Antibody, FITC conjugated

Product Code	CSB-PA009068LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09958
Immunogen	Recombinant Human Furin protein (520-701AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Furin (EC 3.4.21.75) (Dibasic-processing enzyme) (Paired basic amino acid residue-cleaving enzyme) (PACE), FURIN, FUR PACE PCSK3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FURIN