



GABRB3 Antibody, Biotin conjugated

Product Code	CSB-PA009148LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28472
Immunogen	Recombinant Human Gamma-aminobutyric acid receptor subunit beta-3 protein (328-450AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-aminobutyric acid receptor subunit beta-3 (GABA(A) receptor subunit beta-3), GABRB3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GABRB3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.