




GAPDH Antibody

Product Code	CSB-PA009232HA01DXU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17244
Immunogen	Recombinant Cricetulus griseus Glyceraldehyde-3-phosphate dehydrogenase protein (2-333AA)
Raised In	Rabbit
Species Reactivity	Cricetulus griseus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glyceraldehyde-3-phosphate dehydrogenase (GAPDH) (EC 1.2.1.12) (Peptidyl-cysteine S-nitrosylase GAPDH) (EC 2.6.99.-), GAPDH, GAPD
Immunogen Species	Cricetulus griseus
Research Area	Others
Target Names	GAPDH

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: GAPDH antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 36 kDa

observed band size: 36 kDa