



# GAPDH Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA009232HD01DXU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P17244
<b>Immunogen</b>	Recombinant Cricetulus griseus Glyceraldehyde-3-phosphate dehydrogenase protein (2-333AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Cricetulus griseus
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Glyceraldehyde-3-phosphate dehydrogenase (GAPDH) (EC 1.2.1.12) (Peptidyl-cysteine S-nitrosylase GAPDH) (EC 2.6.99.-), GAPDH, GAPD
<b>Immunogen Species</b>	Cricetulus griseus (Chinese hamster) (Cricetulus barabensis griseus)
<b>Research Area</b>	Others
<b>Target Names</b>	GAPDH