



# GCM2 Antibody, HRP conjugated

Product Code	CSB-PA009324LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75603
Immunogen	Recombinant Human Chorion-specific transcription factor GCMb protein (156-451AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chorion-specific transcription factor GCMb (hGCMb) (GCM motif protein 2) (Glial cells missing homolog 2), GCM2, GCMB
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	GCM2